



# VENTURER


888IX

## FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Headlights
- Rear lights
- Indicators
- Hazard lights
- Large battery capacity
- Battery power indicator
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

## SPECIFICATIONS

Overall: 1250mm (L) x 600mm (W) x 1160mm (H)  
 Turning Circle: 1440mm  
 Speed: 9.5km/h  
 Maximum range: Up to 40km  
 Safe climb angle: 9°  
 Maximum climb angle: 9°  
 Motor: 1.0Hp high torque motor (2 Pole)  
 Batteries: 2 x 50 Amp minimum (not supplied)  
 Charger: 8 Amp  
 Tyres (Pneumatic): Ø260mm - Front + Rear  
 Product Weight: 99kg - with batteries  
 Maximum User Weight: 136kg - on flat ground



0° - 136kg  
 6° - 113kg  
 9° - 91kg