

BURGUNDY RED

888WNULLHD

FEATURES

- High back contoured seat with headrest
- Height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Fold up adjustable height / depth footrest
- Mid wheel drive
- Finger control operation
- High / Low speed settings
- Battery power indicator
- Horn
- Ideal for indoor and outdoor use on even ground
- Tie down points for transport

- Attachment options available:
 Swing away legrests
 Elevating swing away legrests
 Stump support

SPECIFICATIONS

Overall: 1060mm (L) x 640mm (W) x 1180mm (H)

Turning Circle: 1090mm

Speed: 6km/h

Maximum range: Up to 35km

Safe climb angle: 8°

Maximum climb angle: 9°

Motor: 2 x 1.34Hp high torque motors (2 Pole)


Batteries: 2 x 75 Amp minimum (supplied)

Charger: 8 Amp

Tyres: Front: Ø180mm / Middle: Ø360mm / Rear: Ø180mm

Product Weight: 126kg - with batteries

Maximum User Weight: 205kg - on flat ground



0° - 205kg
 6° - 150kg
 9° - 113kg