



FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- LCD display
- Indicators / Hazard lights
- Large battery capacity
- Battery power indicator
- USB power outlet
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

- Overall: 1590mm (L) x 730mm (W) x 1285mm (H)
- Turning Circle: 1860mm
- Speed: 10km/h
- Maximum range: Up to 54km
- Safe climb angle: 12°
- Maximum climb angle: 12°
- Motor: 2.0Hp high torque motor (4 Pole)
- Batteries: 2 x 60 Amp minimum (2 x 75 amp supplied)
- Charger: 8 Amp
- Tyres (Pneumatic): Ø330mm - Front + Rear (80/80-8)
- Product Weight: 147kg - with batteries
- Maximum User Weight: 225kg - on flat ground

