SHOPRIDER

G424M

BURGUNDY RED



0° - 136kg 6° - 113kg 9° - 75kg

FEATURES

- High back contoured seat with headrest
- Height / angle adjustable seat
- Lift up & adjustable width / height arm rests
- Fold up adjustable height / depth footrest
- Mid wheel drive
- Finger control operation
- High / Low speed settings
- Battery power indicator
- Horn
- Ideal for indoor and outdoor use on even ground
- Tie down points for transport
- Attachment options available: Swing away legrests Elevating swing away legrests Stump support

SPECIFICATIONS

Overall: 920mm (L) x 620mm (W) x 1120mm (H)

Turning Circle: 950mm

Speed: 8km/h

Maximum range: Up to 35km

Safe climb angle: 8° Maximum climb angle: 9°

Motor: 2 x 0.8Hp high torque motors (2 Pole) Batteries: 2 x 40 Amp minimum (supplied)

Charger: 5 Amp

Tyres: Front: Ø150mm / Middle: Ø260mm / Rear: Ø125mm

Product Weight: 100kg - with batteries

Maximum User Weight: 136kg - on flat ground