

SHOPRIDER®

MOBILITY SCOOTERS & POWERCHAIRS

ROCKY 12

BLUE



REAR
WHEEL DRIVE

FULL LCD
DISPLAY



FEATURES


889WB

- High back contoured seat with headrest
- Wide seat option available
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Emergency hand brake
- Twin headlights (LED)
- Rear lights with brake-lights
- Full LCD display
- Automotive style Indicators / Hazard lights
- Large battery capacity 2 x 100Amps - XT100 Premium
- Battery power indicator
- USB power outlet
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

SPECIFICATIONS

Overall: 1740mm (L) x 750mm (W) x 1310mm (H)
 Std. Seat: 560mm (W) x 540mm (D) - (22" x 21")
 Wide Seat: 820mm (W) x 450mm (D) - (32" x 18")
 Turning Circle: 2705mm
 Speed: 10km/h
 Maximum range: Up to 50km
 Safe climb angle: 9°
 Maximum climb angle: 9°
 Motor: 2Hp high torque motor (4 Pole)
 Batteries: 2 x 12V 100Ah
 Charger: 8 Amp
 Tyres: Ø410mm - Front + Rear
 Product Weight: 186kg
 Maximum User Weight: 350kg - on flat ground

WIDE SEAT
OPTION AVAILABLE



0° - 350kg
 6° - 300kg
 9° - 250kg

