


# RAINRIDER

889XLSBNRR

INJECTION RED

INJECTION BLUE



0° - 160kg  
6° - 110kg  
9° - 90kg  
12° - 70kg

## FEATURES

- Fully enclosed cabin with double sided doors
- Windscreen wipers
- Double pillow comfort high back seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators / Hazard lights
- Large battery capacity
- Battery power indicator
- USB power outlet
- Sliding / Rotating / Height adjustable seat

## SPECIFICATIONS

Overall: 1600mm (L) x 790mm (W) x 1620mm (H)  
 Turning Circle: 1800mm  
 Speed: 10km/h  
 Maximum range: Up to 40km  
 Safe climb angle: 9°  
 Maximum climb angle: 12°  
 Motor: 2.0Hp high torque motor  
 Batteries: 2 x 75 Amp minimum (supplied)  
 Charger: 8 Amp  
 Tyres (Pneumatic): Ø330mm - Front + Rear  
 Product Weight: 195kg - with batteries  
 Maximum User Weight: 160kg - on flat ground