

Shoprider STREAMER

ASK FOR THE BRAND SHOPRIDER



STREAMER

SHOPRIDER[®]
mobility

www.shoprider.com.au

RELIABLE MOBILITY GUARANTEED

Shoprider

STREAMER

SHOPRIDER[®]
www.shoprider.com.au *mobility*

The SHOPRIDER™ **Streamer** is another great powerchair from Shoprider. This budget model offers many of the features found on much more expensive powerchairs. The Streamer comes with a high back adjustable seat, adjustable width & height armrests, flat free tyres and finger control steering. The Streamer is rear wheel drive and combined with two powerful 0.8hp motors it will deliver outstanding performance. This powerchair offers excellent value for money and reliability only a Shoprider can offer.

SPECIFICATIONS	Streamer	FEATURES:
Length	900mm	High back seat
Width	600mm	Sleek body design
Height	1016mm	Smooth operation
Battery Capacity (Minimum Amp)	2 x 36Amp	Variable seat height
#Maximum Range	up to 18km	Adjustable arm rests
Maximum Load Single User*	136Kgs	Controller interchangeable both arms
Charger off board	5Amp	Free wheel device
Motor Size	0.8 Hp	Option of swing away leg rest
Maximum Speed	7.5Kph	Option of elevating leg rest
*Safe Climbing Angle	9 Degrees	Option of Stump Support
*Max Climbing Angle	9 Degrees	
Turning Radius	650mm	
Wheel Size:		
- Front Caster	2 x 150mm	
- Rear drive wheel	2 x 265mm	

Driver weight may exceed weight of unit. Reduce speed when turning

IMPORTANT

Always be sure to read your owner's manual thoroughly before operating your machine.

Driving range varies with the factors shown below

- (1) The weight of the occupant
- (2) Level ground & Ambient temperature
- (3) Battery condition / Type of charger
- (4) Driving Style



Streamer

9' (15%) -91kg
6' (10%) -113kg
0' (0) -136kg

SHOPRIDER[®]
www.shoprider.com.au *mobility*

Your Authorised dealer

For your nearest Distributor call 1800 672 633
Or visit our Web Site at www.shoprider.com.au