



In partnership with  
*Life & Mobility*



METALLIC  
SILVER  
IK8PRS

# PRESTO

## MOBILITY SCOOTER

Part of the **FOR-DE Group Pty Ltd**



# EUROPEAN DESIGN



### FEATURES

- Latest European design features
- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- LED Headlight
- Rear lights
- LCD display
- Indicators / Hazard lights
- Large battery capacity
- Battery power indicator
- USB power outlet
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width armrests
- Removable shopping basket

### WARRANTY

- Structural Frame: **5 years**
- Drive Train Assembly: **18 months**
- Other Parts: **12 months**
- XT 75 Amp Hybrid Gel Batteries: **18 months**
- Victron Premium Charger: **5 years**

**Flip over for Specifications**

#### IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.  
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

P: **08 9248 4180**

E: **orders@for-de.com.au**

ASK FOR THE BRAND

ikonmobility.com.au


IKON Mobility strives to provide the highest quality products possible that are subject to change without notice.




Please contact IKON Mobility if you would like any further information on the Presto Mobility Scooter.

# PRESTO


MOBILITY SCOOTER




**OVERALL LENGTH**  
1340mm




**OVERALL HEIGHT**  
650mm




**OVERALL WIDTH**  
1135mm




**SEAT WIDTH**  
510mm




**SEAT DEPTH**  
480mm




**SEAT HEIGHT**  
450 - 525mm  
(4 adjustments)




**ARMREST LENGTH**  
380mm




**ARMREST WIDTH**  
60mm




**ARMREST HEIGHT  
(FROM SEAT)**  
190mm



**TILLER WIDTH**  
470mm




**TILLER HEIGHT**  
930mm




**TURNING CIRCLE**  
1600mm




**SPEED**  
10km/h




**MAXIMUM RANGE**  
Up to 44km




**SAFE CLIMB ANGLE**  
10°




**MAXIMUM CLIMB ANGLE**  
15°




**MOTOR**  
2.0Hp high torque motor  
(4 Pole)



**BATTERY**  
2 x 75Amp minimum  
(supplied)



**BATTERY CHARGER**  
8 Amp



**TYRES (PNEUMATIC)**  
Ø356mm - Front + Rear



**PRODUCT WEIGHT**  
136kg



**MAXIMUM USER WEIGHT**  
160kg  
(on flat ground)

ASK FOR THE BRAND

ikonmobility.com.au