



IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

IKON Mobility strives to provide the highest quality products possible that are subject to change without notice.



AIRLINE TRAVEL READY



METALLIC GREEN
IK4QPG

METALLIC RED
IK4QPR

TRAVEL
FOLDING MOBILITY SCOOTER

FEATURES

- Folds up for transport and storage
- Height adjustable pull handle
- Safe to take on an aeroplane
- Comfortable contoured seat
- Fold down backrest
- Adjustable height T handlebar
- Finger operated controls
- Battery powered indicator
- Forward / Reverse drive settings
- Variable speed setting
- Rear wheel drive
- Electromagnetic brake
- Rechargeable Lithium-ion battery
- Graduated tru-view Lithium battery gauge

WARRANTY

- Structural Frame: **5 years**
- Drive Train Assembly: **18 months**
- Other Parts: **12 months**
- Lithium Battery: **18 months**
- Charger: **12 months**

Flip over for Specifications

P: **08 9248 4180**

E: **orders@for-de.com.au**


ASK FOR THE BRAND

ikonmobility.com.au




Please contact IKON Mobility if you would like any further information on the eTRAVEL Folding Mobility Scooter.


eTRAVEL FOLDING MOBILITY SCOOTER




OVERALL LENGTH
910mm




OVERALL HEIGHT
870mm




OVERALL WIDTH
420mm




FOLDED LENGTH
710mm




FOLDED WIDTH
420mm




FOLDED HEIGHT
280mm




TURNING CIRCLE
1420mm




SEAT WIDTH
330mm




SEAT DEPTH
250mm



SEAT HEIGHT
390mm
(from footrest)




PULL HANDLE HEIGHT
690 - 1000mm




TILLER HEIGHT
700 - 750mm




SPEED
6.4km/h




MAXIMUM RANGE
Up to 17km




MAXIMUM CLIMB ANGLE
6°




MOTOR
180w high torque motor




BATTERY
24V / 9.9Ah
Lithium-Ion
(supplied)



BATTERY CHARGER
2 Amp
Lithium-Ion
(supplied)



TYRES (FLAT FREE)
Ø165mm - Front
Ø185mm - Rear



PRODUCT WEIGHT
19kg (with battery)
17kg (without battery)



MAXIMUM USER WEIGHT
115kg
(on flat ground)

ASK FOR THE BRAND

ikonmobility.com.au