


NEW PRODUCT RELEASE !



Blue



Grey



Mushroom



Brown



IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

KATIE LEG REST / STOOL

RG63LRBL - BLUE

RG63LRGR - GREY

RG63LRM - MUSHROOM

RG63LRBR - BROWN

FEATURES

Designed to accompany the KATIE Chair range:

- Katie High Back Chair
- Katie Low Back Chair

Can be used as a Leg rest or as a Stool

Available in 4 colours

Padded platform

Height adjustable

Angle adjustable

Strong lightweight construction

Vinyl upholstery

Epoxy powder coated metal framework

Plastic leg tips

SPECIFICATIONS (All dimensions shown are nominal)

Platform Width:

450mm

Platform Depth:

380mm

Platform Height:

325 - 475mm (from floor to top of platform)
25mm adjustment increments

Product Weight:

5kg

Maximum User Weight:

120kg (Evenly distributed)

Please contact **REDGUM Brand** if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090
Part of the **FOR-DE Group** Pty Ltd

redgumbrand.com.au