

FOLDS FOR STORAGE!



PADDED ADJUSTABLE
HEIGHT ARM RESTS

PEARL DELUXE ALUMINIUM WHEELCHAIR

RG18AADRB | RG20AADRB

FEATURES

- Aluminium anti corrosive frame
- Padded adjustable height arm rests
- Adjustable rear handles
- Rear hand brakes fitted on handles
- Adjustable comfort back rest (Velcro)
- Four stage seat adjustment
- Height and angle adjustable foot plates
- Easy to operate lever swing away leg rests
- Manual brake on rear wheels
- Quick disconnect rear wheels
- Anti tip wheels fitted as standard

SPECIFICATIONS

All dimensions shown are nominal and at lowest setting

RG18AADRB

Overall: **670mm (W) x 960mm (H) x 1110mm (L)**

Folded: **400mm (W) x 960mm (H) x 810mm (L)**

Seat: **470mm (W) x 390mm - 430mm (D)**

Armrest Height (from seat): **230mm - 270mm**

Backrest Height (from seat): **410mm**

Footrest Height (from seat): **470mm - 510mm**

Wheels: **Front - 180mm/7" Back - 600mm/24"**

Product Weight: **20kg**

Max. User Weight: **120kg**

RG20AADRB

Overall: **670mm (W) x 960mm (H) x 1110mm (L)**

Folded: **400mm (W) x 960mm (H) x 810mm (L)**

Seat: **510mm (W) x 390mm - 430mm (D)**

Armrest Height (from seat): **230mm - 270mm**

Backrest Height (from seat): **410mm**

Footrest Height (from seat): **470mm - 510mm**

Wheels: **Front - 180mm/7" Back - 600mm/24"**

Product Weight: **22kg**

Max. User Weight: **130kg**



EASY TO OPERATE LEVER
SWING AWAY LEG RESTS

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

Please contact **REDGUM Brand** if you require further information on this product.

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand: 1 Business way Malaga Western Australia 6090
Part of the For-de Group Pty Ltd

redgumbrand.com.au