

**WIRE BASKET INCLUDED!**



**HEIGHT  
ADJUSTABLE  
FOREARM  
GUTTER  
ARMRESTS**



**FLAT FREE  
WHEELS**

**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

# WALKING TUTOR

## RG4310AA

### FEATURES

- Designed to tutor the user into an upright gait increasing balance and stability.
- Height adjustable forearm gutter armrests
- Folds easily for storage and transportation
- Loop style lockable hand brakes
- Flat free wheels
- Includes wire basket and tray

### SPECIFICATIONS

All dimensions shown are nominal and at lowest setting

Overall: (at lowest setting)

**590mm (W) x 700mm (D) x 975 - 1150mm (H)**

Overall Folded: (laid on ground)

Width: **590mm**

Depth: **975mm - 1150mm**

Height: **410mm**

Armrests: **975mm - 1150mm** (from ground)

Armrests: **380mm (W) x 220mm (D)**

Wheels: **ø200mm / 8"**

Product Weight: **14.5kg**

Max. User Weight: **120kg**

Please contact  
**REDGUM Brand** if  
you require further  
information on  
this product.

**FOLDS FOR STORAGE**



**P: 08 9248 4180**

**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

REDGUM brand: 1 Business way Malaga Western Australia 6090  
Part of the **For-de Group Pty Ltd**

**[redgumbrand.com.au](http://redgumbrand.com.au)**