

**BARIATRIC 200KG SWL!**



# TITAN HD WHEELCHAIR

RG22HD



**FOLDS FOR  
TRANSPORT/  
STORAGE**



## FEATURES

- Powder coated steel folding frame
- Padded swing up armrests
- Removable swing away leg rests with adjustable height
- Attendant and occupant hand brakes
- Self propelled quick release wheels
- Depth adjustable seat
- Large swing up foot rest
- Double cross bar seat support
- Foot plate joiner and leg strap included
- Optional elevating leg rests  
(L: **RGWCP1/R: RGWCP2**)
- Optional stump support  
(L: **RGWCP3/R: RGWCP4**)

## SPECIFICATIONS

(All dimensions shown are nominal)

- Overall: W: 810mm x D: 1140mm x H: 955mm
- Folded without leg rests and wheels:  
W: 335mm x D: 900mm x H: 935mm
- Seat: W: 580mm x D: 390mm-435mm x  
H: 480mm
- Backrest: W: 590mm x H: 430mm
- Wheels: (Ø600mm / Ø24" x 1- 3/8")
- Product weight with wheels and leg rests: 20kg
- Wheel Weight: 2.5kg each
- Leg Rests Weight: 0.90kg each
- Max User Weight: 200kg

Please contact **REDGUM Brand** if you require further information on this product



**NOTE:**  
Product comes with attachable foot plate joiner and leg strap.

**P: 08 9248 4180**

**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

1 Business Way Malaga Western Australia 6090  
Part of the **For-de Group Pty Ltd**

**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

**ASK FOR THE BRAND**

**redgumbrand.com.au**