

**BARIATRIC  
WHEELCHAIR!**



Lift up armrests  
for access



Folding frame for  
storage and transport

Removable  
swing away  
leg rests

**IMPORTANT NOTE:**

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

# TRAKKA WHEELCHAIR

RG18TW | RG20TW

## FEATURES

- Bariatric wheelchair - 160kg
- Two sizes available (18" & 20" seat width)
- Folding frame for storage & transport
- Swing up padded armrests
- 45mm thick tyres
- Adjustable seat depth
- Removable swing away leg rests with adjustable height and heel straps
- Padded breathable nylon upholstery
- Manual brake on rear wheels
- Flat free wheels

## SPECIFICATIONS

All dimensions shown are nominal and at lowest setting

### RG18TW

Overall: 650mm (W) x 1070mm (L) x 880mm (H)

Folded: 340mm (W) 780mm x (L) x 880mm (H)

Seat: 450mm (W) x 415mm - 455mm (D) x 470mm (H)

Backrest Height: 410mm - from seat

Footrest Height: 400 - 540mm - from seat

Wheels: Front -  $\phi$  180mm / 7" Back -  $\phi$  570mm / 22"

Product Weight: 20kg

Max. User Weight: 160kg

### RG20TW

Overall: 650mm (W) x 1070mm (L) x 880mm (H)

Folded: 340mm (W) 780mm x (L) x 880mm (H)

Seat: 510mm (W) x 415mm - 455mm (D) x 470mm (H)

Backrest Height: 410mm - from seat

Footrest Height: 400 - 540mm - from seat

Wheels: Front -  $\phi$  180mm / 7" Back -  $\phi$  570mm / 22"

Product Weight: 22kg

Max. User Weight: 160kg

Please contact REDGUM Brand if you require further information on this product.

**P: 08 9248 4180**

**E: [orders@for-de.com.au](mailto:orders@for-de.com.au)**

REDGUM brand: 1 Business way Malaga Western Australia 6090  
Part of the For-de Group Pty Ltd

**[redgumbrand.com.au](http://redgumbrand.com.au)**