



**IDEAL FOR MOVING
SMALL ITEMS
AROUND THE HOUSE**



REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

IMPORTANT NOTE:
It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

TRAY WALKER


RG4520

FEATURES

- Walking frame and transport trolley in one.
- Assist with moving small items around the house. (Food from kitchen to the dining table)
- Ideal for use in confined spaces.
- Height adjustable.
- Removable plastic tray and tub.
- Various tray / tub position settings.
- Twin wheel grey rubber castors.

SPECIFICATIONS

- Overall: **500mm L x 490mm W x 820mm H** (lowest setting)
- Handle Height: **820 - 950mm** (6 settings)
- Tray (Internal): **360mm W x 250mm D x 20mm H**
- Tub (Internal): **330mm W x 230mm D x 100mm H**
- Wheel Diameter: **Ø75mm / 3"**
- Product Weight: **8kg**
- Max. Support Weight: **150kg**
- Max. Tray Load Weight: **10kg**
- Max. Tub Load Weight: **10kg**



Please contact **REDGUM Brand** if you require further information on this product

P: 08 9248 4180

E: orders@for-de.com.au

A: 1 Business way Malaga Western Australia 6090

Part of the **For-de Group Pty Ltd**