

TILT TOP WITH
7 ANGLE ADJUSTMENTS

OVER BED CHAIR TABLE

RG625



HEIGHT
ADJUSTABLE
PLATFORM

WIDTH ADJUSTABLE

OVER BED CHAIR TABLE
CAN BE USED WITH:

- BEDS
- KATIE CHAIRS
- LIFT CHAIRS



Blue High Back Katie Chair
(RG630BL) not included*



FEATURES

- Large working surface
- Height adjustable split platform
- Width adjustable base
- Tilt top with 7 angle adjustments
- Fixed side platform
- Twin wheel locking castors
- Epoxy powder coated metal frame
- Chrome mild steel upright
- Can be flat packed for storage and transport

SPECIFICATIONS

(All dimensions shown are nominal)

- Large Platform: **750mm (L) x 400mm (D)**
- Side Platform: **150mm (L) x 400mm (D)**
- Platform Height: **700 - 980mm**
- Base Frame: **540 - 765mm (W - internal) x 100mm (H)**
- Castors: **Ø50mm plastic twin wheel locking castors**
- Product Weight: **8kg**
- Maximum Weight Loading: **10kg - evenly distributed**

Please contact **REDGUM Brand** if you require further information on this product.

P: **08 9248 4180**
E: **orders@for-de.com.au**

ASK FOR THE BRAND

redgumbrand.com.au