



IDEAL FOR INCONTINENCE MATTRESS PROTECTION

FLUID INGRESS 150ML P/H PENETRATION RATE

CAN BE POSITIONED HORIZONTALLY & VERTICALLY FOR SPECIFIC COVERAGE



Top Covering:
100% Polyester Microfibre
Underside:
Non woven 100% Polypropylene

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

MATTRESS PROTECTIVE COVER

RG66MC

FEATURES

- Ideal for use on any size mattress
- Protects your mattress from accidental stains
- Spill barrier between mattress and user
- Ideal for incontinence purposes*
- Clean and hygienic - washable
- Can be placed width ways to tuck under mattress or length ways to offer more coverage
- High quality fabrics / materials

*Fluid ingress: 150ml p/h penetration rate

SPECIFICATIONS

- Overall: W: 1980mm x D: 895mm x H: 1mm
- Product weight: 0.58kg

CLEANING

- Spot clean any spills as soon as possible using a non-abrasive cleaner or mild detergent with a soft cloth.
- Use a soft bristled brush on hard stains.
- Wash separately (avoid cottons)
- Warm (40°C or less) or cool wash with detergent.
- Hang to dry or tumble dry on low heat.
- Do not use fabric softeners or bleach.
- Care instructions can also be found on the label sewn onto the underside of the cover.

Please contact **REDGUM Brand** if you require further information on this product

P: 08 9248 4180

E: orders@for-de.com.au

1 Business Way Malaga Western Australia 6090
Part of the **For-de Group Pty Ltd**