

BURGUNDY RED



QT8-3

FEATURES

- Scooter pulls apart for storage and transport
- Comfortable mid back contoured seat
- Fold down adjustable steering tiller
- Straight handle with finger operated controls
- Rear wheel drive
- Electromagnetic brake
- Battery power indicator
- Horn
- Removable rotating / height adjustable seat
- Lift up & adjustable width arm rests
- Lifting handles makes picking up the scooter easy
- Removable shopping basket

SPECIFICATIONS

Overall: 960mm (L) x 560mm (W) x 985mm (H)

Turning Circle: 830mm

Speed: 6km/h

Maximum range: Up to 18km

Safe climb angle: 6°

Maximum climb angle: 6°

Motor: 0.3Hp high torque motor (2 Pole)

Batteries: 2 x 12 Amp minimum (supplied)

Charger: 2 Amp

Tyres (Flat free): Ø200mm - Front + Rear

Maximum User Weight: 113kg - on flat ground

Product Weight: 41kg - with batteries

Battery Weight: 5.6kg each

Seat Weight: 8.3kg

Rear Transaxle Weight: 9.5kg

Main Chassis Weight: 12kg

